

LOTUS CATHETER STABILIZER INSTRUCTIONS FOR USE

When the indications for the Lotus Catheter use as a Foley Balloon Catheter or as a Maelcot Catheter. That is indwelling for up to 30 days. It would be preferable to use the Lotus Catheter Stabilizer instead of using the standard adhesive tapes to stabilize the bellows to the inner side of the patient's thigh. . The Lotus Catheter adhesive pad is hypo allergenic. It also has a very strong bonding with the skin.

!- After identifying exactly the position of the Lotus Catheter inside the urinary bladder. The Lotus catheter site on the inner side of the thigh is cleaned, prepped, and dried. In male patients shaving the site is important as the adhesive used has strong bondage.

2- The adhesive part of the Lotus Catheter Stabilizer may be tailored in size to fit a pediatric or a young small patient.

3- The Lotus catheter bellows is placed firmly in the bed of the Stabilizer.

4- The moving covering lid has a latch. By closing this lid, the latch would firmly lock into the recipient of the latch on the inner side of the stabilizer bed.

5- To remove the Lotus Catheter Stabilizer. Press the covering moving lid firmly towards the bellows to disengage the latch from the latch recipient. Then release the bellows from the stabilizer's bed

6- In difficult cases to open the locked covering lid. Insinuate a small, closed nursing scissor between the moving covering lid and the bed of the Stabilizer. Twist the closed scissors five degrees or so clockwise. This will unlock the latch from the latch recipient. Then remove the bellows from the bed of the Lotus Stabilizer.

7- In patients with sensitive skin. It would be advantageous to use alcohol wipes slowly and patiently to totally remove the adhesive patch of the Lotus Stabilizer with minimal or no discomfort to the patient.